

Report no.	Product	Test Standard	Category	Concentration	Contact Time	Notes
E166 *	Ultraseptin classic	EN13697:2015	Bacteria	0.5%	5 minutes	*Replaced E123 to include successful Fungi testing – other results are identical to E123
			Yeast	0.5%	5 minutes	
			Fungi	2.0%	60 minutes	
E124	Ultraseptin classic	EN13727:2015	Bacteria	0.5%	5 minutes	-
E165*	Ultraseptin classic	EN13624:2013	Yeast	0.5%	5 minutes	*Replaced E125 to include successful Fungi testing – other results are identical to E125
			Fungi	2.0%	60 minutes	
E126	Ultraseptin classic	EN14561:2006	Bacteria	1.0%	5 minutes	-
E127	Ultraseptin classic	EN14562:2006	Yeast	0.5%	5 minutes	-
			Fungi	2.0%	30 minutes	
E128	Ultraseptin classic	EN14476:2013+A1:2015	Virus (Polio)	1.0%	30 minutes	-
				2.0%	15 minutes	
				2.0%	15 minutes	
E131	Ultraseptin (All)	EN14476:2013+A1:2015	Virus (Polio)	1.0%	30 minutes	Cornerstone efficacy testing - for information only. Related to E128-E130 to allow report transcription
E132	Ultraseptin classic	EN14348:2005	Mycobacteria	1.0%	30 minutes	-
				2.0%	15 minutes	
E133	Ultraseptin classic	EN14563:2009	Mycobacteria	0.5%	5 minutes	-
				1.0%	5 minutes	
				2.0%	30 minutes	

Report no.	Product	Test Standard	Category	Concentration	Contact Time	Notes
E134	Ultraseptin classic	EN14476:2013+A1:2015	Virus (MNV)	0.50%	2 minutes	-
E135	Ultraseptin classic	EN14476:2013+A1:2015	Virus (Adeno)	0.50%	5 minutes	-
E136	Ultraseptin classic	EN14476:2013+A1:2015	Virus (Virucidal cornerstone efficacy Polio)	1.0%	30 minutes	-
				2.0%	15 minutes	
E139	Ultraseptin (All)	EN14476:2013+A1:2015	Virus (Virucidal cornerstone efficacy Adeno)	1%	30 minutes	Cornerstone efficacy testing - for information only. Related to E135, E137, E138 to allow report transcription
E142	Ultraseptin (All)	EN14476:2013+A1:2015	Virus (Virucidal cornerstone efficacy MNV)	1%	30 minutes	Cornerstone efficacy testing - for information only. Related to E134, E140, E141 to allow report transcription